


# Easy Enrichment of VOC from Papers Before & After Printing with MonoTrap



MonoTrap is a hybrid novel adsorbent that combines large surface area and properties of silica gel, activated carbon, and ODS. Utilizing the surface of porous silica, it can be used as a tool that makes use of the large surface area and the adsorbent effect caused by the inclusion of activated carbon for high collection efficiency and for short time and sensitive analysis. MonoTrap DCC18 (with activated carbon) was used to analyze VOCs in new vehicles. Activated carbon is effective in collecting polar and low boiling components. It is also easy to collect by simply placing a small collection material inside the car room. Capture concentration can also be used as usual, and simple and stress-free enrichment analysis can be performed.

## Protocol



### GC Conditions

<b>System</b>	: GC/MS
<b>Column</b>	: InertCap Pure-WAX 0.25 mm I.D. × 30 m df = 0.25 µm
<b>Col. Temp.</b>	: 40 °C (5 min) - 4 °C /min - 250 °C
<b>Carrier Gas</b>	: He 1 mL/min
<b>Injection</b>	: Splitless 250 °C
<b>Detection</b>	: MS Scan ( <i>m/z</i> 35-500)
<b>Sample Size</b>	: 1 µL



GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

**GL Sciences Inc. Japan**  
22-1 Nishishinjuku 6-chome  
Shinjuku-ku, Tokyo  
163-1130, Japan

Phone: +81-3-5323-6620  
Fax: +81-3-5323-6621  
Email: [world@glsc.co.jp](mailto:world@glsc.co.jp)  
Web: [www.glsciences.com](http://www.glsciences.com)

**GL Sciences Inc. USA**  
4733 Torrance Blvd. Suite 255  
Torrance, CA 90503  
USA

Phone: +1-310-265-4424  
Fax: +1-310-265-4425  
Email: [info@glsciencesinc.com](mailto:info@glsciencesinc.com)  
Web: [www.glsciencesinc.com](http://www.glsciencesinc.com)

**GL Sciences B.V.**  
Dillenburgstraat 7C  
5652AM, Eindhoven  
The Netherlands

Phone: +31-40-254-95.eu  
Web: [www.glsciences.eu](http://www.glsciences.eu)

**GL Sciences (Shanghai) Limited**  
Tower B, Room 2003  
Far East International Plaza  
No.317 Xianxia Road, Changning District  
Shanghai, China 200051

Phone: +86-21-62782272  
Email: [contact@glsciences.com.cn](mailto:contact@glsciences.com.cn)  
Web: [www.glsciences.com.cn](http://www.glsciences.com.cn)



**International Distributors**

Visit our Website at [www.glsciences.com/distributors](http://www.glsciences.com/distributors)