


Analysis of Residual Solvent Components in Pharmaceutical Products – Using InertCap Pure-WAX

This application is an example of the analysis of components subject to the Guideline for Residual Solvents in Drugs.

Depending on the level of toxicity, this guideline classifies solvents as Class 1, Class 2, Class 3, and those for which no appropriate toxicity data are available. For these solvents, simultaneous analysis of the components was performed using GC/FID.

Example: measurement of mixed sample



Sample name and retention index

InertCap Pure WAX					
Component	RI	Component	RI	Component	RI
1. <i>n</i> -Pentane	500	25. 1,2-Dimethoxyethane	931	48. <i>p</i> -Xylene	1143
2. <i>n</i> -Hexane	600	26. 2-Propanol(Isopropyl alcohol)	931	49. 1-Butanol	1151
3. Diethyl ether	632	27. Dichloromethane	937	50. <i>m</i> -Xylene	1151
4. 2,2,4-Trimethylpentane	670	28. Ethanol	937	51. Nitromethane	1160
5. Diisopropyl ether	670	29. Methyl isopropyl ketone	937	52. 2-Methoxyethanol	1186
6. <i>tert</i> -Butyl methyl ether	691	30. Benzene	948	(Methyl cellosolve)	
7. <i>n</i> -Heptane	700	31. Propionaldehyde diethyl acet	963	53. Cumene	1186
8. 1,1-Dimethoxymethane	729	32. <i>n</i> -Propyl acetate	984	54. Pyridine	1191
9. Cyclohexane	729	33. <i>cis</i> -1,2-Dichloroethylene	1001	55. <i>o</i> -Xylene	1195
10. 1,1-Dichloroethylene	729	34. Trichloroethylene	1004	56. 3-Methyl-1-butanol	1217
11. Methylcyclohexane	777	35. Acetonitrile	1012	(Isoamyl alcohol)	
12. Acetone	819	36. 4-Methyl-2-pentanone(MIBK)	1016	57. Chlorobenzene	1225
13. 2,2-Dimethoxypropane	830	37. Isobutyl acetate	1022	58. 2-Ethoxyethanol	1231
14. Methyl acetate	830	38. Chloroform	1031	59. 1-Pentanol(Amyl alcohol)	1259
15. Ethyl formate	830	39. 2-Butanol	1031	60. N,N-Dimethylformamide	1335
16. Tetrahydrofuran	866	40. 1-Propanol	1044	61. Anisole	1354
17. <i>trans</i> -1,2-Dichloroethylene	866	41. Toluene	1049	62. N,N-Dimethylacetamide	1413
18. 2-Methyltetrahydrofuran	876	42. 1,4-Dioxane	1069	63. Acetic acid	1448
19. 1,1,1-Trichloroethane	891	43. 1,2-Dichloroethane	1074	64. 1,2,3,4-Tetrahydronaphthalene	1541
20. Carbon tetrachloride	891	44. <i>n</i> -Butyl acetate	1083	65. Dimethyl sulfoxide(DMSO)	1574
21. Ethyl acetate	897	45. 2-Hexanone(MBK)	1091	66. Ethylene glycol	1631
22. Methanol	903	46. 2-Methyl-1-propanol	1100	67. N-methyl-2-pyrrolidone	1681
23. 2-Butanone(MEK)	908	(Isobutyl alcohol)		68. Formamide	1787

Xylene used is a mixture of *m*-Xylene, *p*-Xylene, *o*-Xylene, and Ethylbenzene.

Retention indices are...

This value is based on the number of carbons in the straight-chain alkanes and is calculated using the retention time of each component.

In this application note, a temperature-rise analysis was made. The formula is shown below.

$$\text{Retention index} = 100 \times \frac{\text{TR} - \text{tR}(Z)}{\text{TR}(Z+1) - \text{tR}(Z)} + 100 \times Z$$

TR = retention time of the target compound
 T R(Z) = retention time of straight-chain alkanes that precede the components of interest
 TR(Z+1) = retention time of straight-chain alkanes emerging after the component of interest.
 Z = number of carbons in the straight-chain alkane with hydrocarbon retention time tR(Z)

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome
 Shinjuku-ku, Tokyo
 163-1130, Japan

Phone: +81-3-5323-6620
 Fax: +81-3-5323-6621
 Email: world@glsc.co.jp
 Web: www.glsciences.com

GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255
 Torrance, CA 90503
 USA

Phone: +1-310-265-4424
 Fax: +1-310-265-4425
 Email: info@glsciencesinc.com
 Web: www.glsciencesinc.com

GL Sciences B.V.

Dillenburgstraat 7C
 5652AM, Eindhoven
 The Netherlands

Phone: +31-40-254-9531
 Email: info@glsciences.eu
 Web: www.glsciences.eu

GL Sciences (Shanghai) Limited

Tower B, Room 2003
 Far East International Plaza
 No.317 Xianxia Road, Changning District
 Shanghai, China 200051

Phone: +86-21-62782272
 Email: contact@glsciences.com.cn
 Web: www.glsciences.com.cn



International Distributors

Visit our Website at www.glsciences.com/distributors