


Analysis of Bufosteroids in Senso In Accordance with the Japanese Pharmacopoeia

Bufotoxin refers to crude drugs collected from the venomous glands of a toad (e.g., *Bufo bufo garagarizans* Cantor).

The Japanese Pharmacopoeia HPLC method to quantify bufo steroids in senso employs an ODS-column. The following system suitability items are specified.


(T. Fukaya)

Structural formula



Structures are created using Chemistry 4-D Draw which is provided by ChemInnovation Software, Inc.

Example: Measurement of standard



System suitability test

When analyzed under the above HPLC conditions,

1. Bufalin, cinobufagin, resibufogenin and indomethacin (internal standard) are eluted in this order.
2. Each peak is completely separated.

Our results

Elution order : Satisfy the provisions
Resolution : Peak 1 - 2 : 11.8
Peak 2 - 3 : 14.6
Peak 3 - 4 : 17.8



Complete separation






Figure 1: Calibration curve of bufosteroids.

*Refer to Technical Note No. 38 for the resolutions.

Analysis of Drugs Containing Commercially Available Bufotoxin

Example of pretreatment



GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome
Shinjuku-ku, Tokyo
163-1130, Japan

Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
USA

Phone: +1-310-265-4424
Fax: +1-310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

GL Sciences B.V.

Dillenburgstraat 7C
5652AM, Eindhoven
The Netherlands

Phone: +31-40-254-9531
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (Shanghai) Limited

Tower B, Room 2003
Far East International Plaza
No.317 Xianxia Road, Changning District
Shanghai, China 200051

Phone: +86-21-62782272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn



International Distributors

Visit our Website at www.glsciences.com/distributors