

Scopolia extract powder is derived from *Scopolia japonica* Maxi. It is a crude drug prepared by drying the roots and extracting the active ingredients. It contains approximately 0.1 % scopolamine and hyoscyamine in the dried product (atropine in this racemic mixture). It is used in combination with over-the-counter stomach drugs as an antispasmodic and analgesic because it reduces the secretion of stomach juice and excessive motility of the gastrointestinal tract, and prevents the transmission of pain.

The Japanese Pharmacopoeia employs an HPLC method with ODS-column for the determination of scopolamine and hyoscyamine in scopolamine extract powder. The following system suitability items are specified.

(T. Fukaya)

Example: Measurement of standard

HPLC conditions

Column : Inertsil ODS-SP
(5 μm , 150 x 4.6 mm I.D.)

Eluent : A) CH_3CN
B) Phosphate buffer
A/B = 10/90, v/v


Flow rate : 0.9 mL/min.

Column temperature : 30 $^\circ\text{C}$

Detected : UV 210 nm

Injection volume : 10 μL

* phosphate buffer: Dissolve a 6.8 g of potassium dihydrogen phosphate in water 900 mL, add triethylamine 10 mL to adjust the pH 3.5 with phosphoric acid, and add water to make a 1000 mL.



System suitability test

When analyzed under the above HPLC conditions,

1. Scopolamine, atropine, and the internal standard are eluted in this order.
2. The resolution of scopolamine and atropine is 11 or greater
3. The resolution between atropine and the internal standard is 4 or greater.

Our results

Separation order: meets specifications

Resolution: scopolamine-atropine: 16.8

Atropine-internal standard: 7.1






Figure 1 Calibration curve for scopolamine and hyoscyamine.

Analysis of scopolamine and hyoscyamine in JP Scopolia Extract

Example of pretreatment




Overview of Resolution

Resolution is a measure of the separation between two peaks and is derived from the retention time and width of two peaks. The greater the number, the greater separation of the two peaks. To calculate resolution, the following formula is used:


$$R_s = 1.18 \times \frac{(t_{R2} - t_{R1})}{(W_{0.5h1} + W_{0.5h2})}$$

t_{R1} , t_{R2} : Retention time of components 1, 2 ($t_{R2} > t_{R1}$)
 $W_{0.5h1}$, $W_{0.5h2}$: peak half-width of components 1, 2
 (Half-width: peak width at peak half height)



Generally, it is considered that determination is possible with a resolution of 1 or more, and complete separation is considered to have been achieved with a resolution of 1.5 or more.


Resolution tends to increase with increased distance in the retention times of the two peaks and with narrower peak widths.




The degree of separation is influenced by the length and type of column, composition of the eluent, column temperature, etc., but is particularly influenced by the type of column. As an example, the results of analyzing the same samples on three ODS columns with different binding modes and different binding amounts of the ODS groups are shown below.

The separation behavior of the same "ODS" column is significantly different, indicating that column selection is extremely important.


Commercial column A



Commercial column B



Inertsil ODS-SP



GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome
Shinjuku-ku, Tokyo
163-1130, Japan

Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
USA

Phone: +1-310-265-4424
Fax: +1-310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

GL Sciences B.V.

Dillenburgstraat 7C
5652AM, Eindhoven
The Netherlands

Phone: +31-40-254-9531
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (Shanghai) Limited

Tower B, Room 2003
Far East International Plaza
No.317 Xianxia Road, Changning District
Shanghai, China 200051

Phone: +86-21-62782272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn



International Distributors

Visit our Website at www.glsciences.com/distributors