


This note describes a determination method for non-volatile putrefactive amines in food using pre-column derivatization.

Putrefactive amines in food are produced by microorganisms. Many of the biogenic amines are known to be responsible for food-poisoning. To determine these amines, HPLC systems coupled with pre-column or post-column derivatization methods are frequently used.

Based on the inspection guideline by Japan Food Hygiene Association, pre-column derivatization method using DNS-Cl (5-dimethylaminonaphthalene-1-sulfonyl chloride) was adopted in this note. The putrefactive amines were successfully fluorescently derivatized, and the calibration curves showed excellent linearity. (K.Suzuki)


Chemical Structures



Structures are created using Chemistry 4-D Draw which is provided by ChemInnovation Software, Inc.


A Chromatogram Obtained from Standard Solution

Standard solution derivatized by DNS-Cl was analyzed.



Conditions

Column	Inertsil ODS-SP (5 μ m, 250 x 4.6 mm I.D.) Cat.No. 5020-02746
Guard column	Cartridge Guard Column E ODS-SP (5 μ m, 10 x 4.0 mm I.D.) Cat.No. 5020-08520
Eluent	A) CH ₃ CN B) H ₂ O A/B= 65/35, v/v, 1.0 mL/min
Column Temp.	40 °C
Detection	FL Ex 325 nm, Em 525 nm
Injection Vol.	10 μ L




The calibration curves for the putrefactive amines

Examples of Sample Pretreatment


Fish samples

Soy sauce samples




*Before loading the sample on SPE cartridge, conditioning was carried out with 10 mL of methanol followed by 10 mL of water.

Extract from fish




Extract from soy sauce



Derivatization reaction by DNS-Cl

DNS-Cl reacts with primary and secondary amines. As a result, fluorescent derivatives are generated.



GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome
Shinjuku-ku, Tokyo
163-1130, Japan

Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
USA

Phone: +1-310-265-4424
Fax: +1-310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

GL Sciences B.V.

Dillenburgstraat 7C
5652AM, Eindhoven
The Netherlands

Phone: +31-40-254-9531
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (Shanghai) Limited

Tower B, Room 2003
Far East International Plaza
No.317 Xianxia Road, Changning District
Shanghai, China 200051

Phone: +86-21-62782272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn



International Distributors

Visit our Website at www.glsciences.com/distributors