

Separation modes such as anion-exchange, ion-exclusion, and reversed-phase have been used in the analysis of organic acids by HPLC. In recent years, more and more cases have been analyzed using ion-exclusion modes with columns dedicated to organic acids.

However, in order to improve the separation in ion exclusion mode, it is necessary to lengthen the analysis time.


In this study, an ion exchange column was used. Analyses were conducted using a combination of "ion exclusion mode" and "hydrophobic interaction mode" with a reversed-phase column. The use of two general-purpose-size columns with different separation modes (250 x 4.6 mm I.D.) enables the rapid simultaneous analysis of 8 components including phosphoric and organic acids.

### Comparison with the previous report (Technical Note No. 1)

The data in Technical Note No. 1 showed that the separation of citric acid, tartaric acid, lactic acid, and formic acid was difficult, but improvement was made.

### \*\* Benefits \*\*

1. Combining the ion exclusion mode with reversed-phase mode allows the separation of 8 organic acids!
2. Postlabelling with BTB (bromothymol blue) reagent eliminates contaminants!
3. Two columns of general purpose size (4.6 x 250 mm I.D.) are used for the separation, which is inexpensive and offers fast to analysis!



### HPLC conditions


Mobile phase : 3mM HClO<sub>4</sub>  
 Reaction solvent : 0.1mM BTB + 30mM Na<sub>2</sub>HPO<sub>4</sub>  
 Columns : Inertsil Ph-3 (250 × 4.6 mm I.D.)  
           + Inertsil CX (250 × 4.6 mm I.D.)  
 Mobile phase flow rate : 0.5 mL/min  
 Reaction solvent flow rate : 0.5 mL/min  
 Temperature : 35 °C  
 Detector : UV-Visible at 440 nm  
 Injection volume : 10 µL

Sample (standard sample)

1. Phosphoric acid (phosphate) (1 mg/mL)
2. Tartaric acid (tartaric acid)(1mg/mL)
3. Malic acid (malate) (1 mg/mL)
4. Formic acid (formic acid) (1 mg/mL)
5. Citric acid (citric acid) (1 mg/mL)
6. Lactic acid (lactic acid) (1 mg/mL)
7. Acetic acid (acetic acid) (1 mg/mL)
8. Succinic acid (succinate) (1 mg/mL)
9. Pyroglutamic acid (pyroglutamic acid) (1 mg/mL)

### \*Detection Lower Limit and Calibration Curve\*

Compound	Concentration (µg/mL)	R <sup>2</sup>
Phosphoric acid	12	0.998
Tartaric acid	16	1.000
Malic acid	13	0.999
Formic acid	10	0.999
Citric acid	15	1.000
Lactic acid	27	1.000
Acetic acid	19	0.995
Succinic acid	15	0.999
Pyroglutamic acid	43	0.999



# Features of organic acid analysis and introduction of analysis system

## \*Features

### To reduce the total cost of organic acid analysis

The use of inexpensive silica columns reduces the cost of the separation columns.

### High resolution

Inertsil CX is a silica gel-based column, the CX media has less diffusion with a greater number of theoretical plates per fixed volume than a polymer-based column.


### Reducing the analysis time

The time required for a single analysis can be reduced to 25 minutes compared to the use of a typical organic acid-column.

Silica-based columns with higher pressure resistance can be used to shorten the analysis time by increasing the flow rate.

## \*Ion Exclusion Mode\*

● Neutral ⊖ weak acid ⊕ Strong acid



The ion exclusion mode uses the strong and weak repulsive force of the ion exchange group on the packing material to separate the target analytes on the basis of ionic strength.

### Determinants of separation time

#### 1. [Penetration into the pores]

The penetration power into the pore is determined by the size of the charge (pKa value).

Higher pKa = greater penetration

#### 2. [Electrostatic exclusion due to the negative charge of the packing material]



Molecules with a large negative charge (small pKa) and are subject to large electrostatic exclusion and elute more quickly.

Small pKa = fast elution

Because molecules are separated by ionic strength, The stronger the anionic strength, the faster the elution.  
Elution order

## \*BTB (bromothymol blue) method\*

Because detection using the BTB method is made using a visible wavelength at 440 nm, samples such as crops, dressings, and soy sauce, which contain many contaminants, can also be analyzed  
See Technical note No. 1 - Highly Selective Assays




Absorption spectrum

The use of BTB reagent results in low absorption at 440 nm for weakly alkaline mobile phase only and a high absorption for weakly acidic samples.

This difference is used to detect the peaks of interest.

# Example of Analysis

Organic acids absorb at UV 210 nm, but when brewed foods such as red wine and beer are analyzed at low wavelengths using UV 210 nm, many contaminants are detected, making it difficult to quantify only the target analyte. However, the BTB method allows selective detection of organic acids because it shifts the detection wavelength to 440 nm.



## Organic Acid Analysis in Red Wine

### Pretreatment

- Sample
  - 1 mL
- Dilution
  - H<sub>2</sub>O 10 mL
- filtration
  - 0.45 μm GL Chromatodisk
- HPLC


1. Phosphoric acid (phosphate)
2. Tartaric acid (tartaric acid)
3. Malic acid (malate)
4. Formic acid (formic acid)
5. Citric acid (citric acid)
6. Lactic acid (lactic acid)
7. Acetic acid (acetic acid)
8. Succinic acid (succinate)
9. Pyroglutamic acid (pyroglutamic acid)

## Organic Acid Analysis in Beer


### Pretreatment

- Sample
  - 1 mL
- Dilution
  - H<sub>2</sub>O 10 mL
- Filtration
  - 0.45 μm GL Chromatodisk
- HPLC

1. Phosphoric acid (phosphate)
2. Tartaric acid (tartaric acid)
3. Malic acid (malate)
4. Formic acid (formic acid)
5. Citric acid (citric acid)
6. Lactic acid (lactic acid)
7. Acetic acid (acetic acid)
8. Succinic acid (succinate)
9. Pyroglutamic acid (pyroglutamic acid)



## Flow diagram



## Analytical column

Inertsil Ph-3 5  $\mu\text{m}$  250 x 4.6 mm I.D.

Cat.No. 5020-01921

Inertsil CX 5  $\mu\text{m}$  250 x 4.6 mm I.D.

Cat.No. 5020-07146

	Product name	Maternal silica gel			Chemical Modification			
		Particle size	Surface area	Pore size	Bonded group	End-capping	Carbon load	Ion Change capacity
Ph	Inertsil Ph-3	5 $\mu\text{m}$	450 $\text{m}^2/\text{g}$	100 $\text{\AA}$	Phenyl group	Non	10 %	-
CX	Inertsil CX	5 $\mu\text{m}$	450 $\text{m}^2/\text{g}$	100 $\text{\AA}$	Benzenesulfonic group	Non	14 %	0.5 meqv/g

Aquatic 25A 0.45  $\mu\text{m}$  filter  
Cat. No. 5040-28512

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

**GL Sciences Inc. Japan**

22-1 Nishishinjuku 6-chome  
Shinjuku-ku, Tokyo  
163-1130, Japan

Phone: +81-3-5323-6620  
Fax: +81-3-5323-6621  
Email: [world@glsc.co.jp](mailto:world@glsc.co.jp)  
Web: [www.glsciences.com](http://www.glsciences.com)

**GL Sciences Inc. USA**

4733 Torrance Blvd. Suite 255  
Torrance, CA 90503  
USA

Phone: +1-310-265-4424  
Fax: +1-310-265-4425  
Email: [info@glsciencesinc.com](mailto:info@glsciencesinc.com)  
Web: [www.glsciencesinc.com](http://www.glsciencesinc.com)

**GL Sciences B.V.**

Dillenburgstraat 7C  
5652AM, Eindhoven  
The Netherlands

Phone: +31-40-254-9531  
Email: [info@glsciences.eu](mailto:info@glsciences.eu)  
Web: [www.glsciences.eu](http://www.glsciences.eu)

**GL Sciences (Shanghai) Limited**

Tower B, Room 2003  
Far East International Plaza  
No.317 Xianxia Road, Changning District  
Shanghai, China 200051

Phone: +86-21-62782272  
Email: [contact@glsciences.com.cn](mailto:contact@glsciences.com.cn)  
Web: [www.glsciences.com.cn](http://www.glsciences.com.cn)

**International Distributors**

Visit our Website at [www.glsciences.com/distributors](http://www.glsciences.com/distributors)