

How to use preparative HPLC - Part 1


Evaluation of loading amount

In preparative HPLC, it is quite important to investigate loading amount per one injection because required times of injection and suitable column dimension can be estimated from loading amount.

In this note, how to evaluate loading amount using analytical column (4.6 mm I.D.) is described.

(K. Kanno)

Conditions of preparative HPLC are generally optimized as the following flowchart.



<What is important to evaluate loading amount>

Main purpose of preparative HPLC is not quantitation but purification of sample. Therefore, good peak shape is not necessarily required. It does not matter, either, if injected amount exceeds detector's upper limit of linearity. Following two points are important to evaluate loading amount;

- 1) All amount of the target compound is well retained on column
- 2) Target compound is successfully separated from other peaks

1) All amount of the target compound is well retained on column


① Prepare sample solution at as high concentration as possible.

② Begin with small amount, and then increase injection volume step-by-step.

In the figure shown below, injection volume was changed as follows; 10 μL \rightarrow 50 μL \rightarrow 100 μL \rightarrow 250 μL . It should be efficient to evaluate peak shape every time increasing injection volume from 2 to 10 times.

③ Peak shape gets worse as injection volume is increased. Injection volume at which peak shape begins to be deteriorated is the maximum amount under the conditions.

In the figure shown below, the peak began to lose its shape when 100 μL of sample solution was injected. Peak was split into two with 250 μL injection. Therefore, it can be said that upper limit of injection volume (amount) under the condition should be around 100 μL (10 mg).



In case the maximum injection volume obtained with 4.6-mm I.D. column is 100 μL as the figure shown above, the maximum injection volume of preparative column is estimated as follows. Detail of scaling up is described in LC Technical Note No.134.

| Inner diameter (I.D.) of column | 4.6 mm | 10 mm | 14 mm | 20 mm | 30 mm |
|---------------------------------|-------------------|-------------------|-------------------|---------|---------|
| Relative cross-sectional area | 1 | 4.7 | 9.3 | 18.9 | 42.5 |
| The maximum injection volume | 100 μL | 470 μL | 930 μL | 1.89 mL | 4.25 mL |


2) Target compound is successfully separated from other peaks

- ① Prepare sample solution at as high concentration as possible.
- ② Begin with small amount, and then increase injection volume step-by-step.

In the figure shown below, injection volume was changed as follows; 50 μL \rightarrow 100 μL \rightarrow 200 μL \rightarrow 400 μL . It should be efficient to evaluate separation of compounds every time increasing injection volume from 2 times to 10 times.

③ The separation gets worse as injection volume is increased. Peak resolution (R_s) is quite useful to evaluate whether sufficient separation is achieved. Target value of R_s varies depending on required purity or how to fractionate the peak, while baseline separation ($R_s > 1.5$) is ideal. R_s value of between 1.2 and 1.5 is often sufficient when fractionation is carried by detecting peak height or peak slope.

In the figure shown below, it can be said that upper limit of injection volume (total loading amount) should be around 200 μL (18 mg) in case the target value of R_s is set at 1.2.



In case the maximum injection volume obtained with 4.6 mm I.D. column is 200 μL as the figure shown above, the maximum injection volume of preparative column is estimated as follows. Detail of scaling up is described in LC Technical Note No.134.

| Inner diameter (I.D.) of column | 4.6 mm | 10 mm | 14 mm | 20 mm | 30 mm |
|---------------------------------|-------------------|-------------------|---------|---------|--------|
| Relative cross-sectional area | 1 | 4.7 | 9.3 | 18.9 | 42.5 |
| The maximum injection volume | 200 μL | 940 μL | 1.86 mL | 3.78 mL | 8.5 mL |

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome
Shinjuku-ku, Tokyo
163-1130, Japan

Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
USA

Phone: +1-310-265-4424
Fax: +1-310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

GL Sciences B.V.

Dillenburgstraat 7C
5652AM, Eindhoven
The Netherlands

Phone: +31-40-254-9531
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (Shanghai) Limited

Tower B, Room 2003
Far East International Plaza
No.317 Xianxia Road, Changning District
Shanghai, China 200051

Phone: +86-21-62782272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn



International Distributors

Visit our Website at www.glsciences.com/distributors