

Detection and determination of 11-nor-9-carboxy-Delta 9-carboxylic acid (THCA), which is a metabolite of principal psychoactive constituent of cannabis plant, is necessary to confirm cannabis use from urine. Sample pretreatment method suitable for handling low volume and low concentration sample has been desired for the determination because THCA

concentration in urine is in ppb level.

In this note, MonoSpin, which provides high recovery even from sample of less than 100  $\mu\text{L}$ , was used for the sample pretreatment, and the purified solution was injected into LC/MS/MS system. The results showed good linearity, recovery, and reproducibility.



(Y. Yui and S. Ota)

### Evaluation of purification by MonoSpin C18-CX (analyzed by HPLC-UV system)

<Before purification>

<After purification>

1. Carboxy-THC (10  $\mu\text{g/mL}$ )



#### Conditions

Column : Inertsil ODS-3 (5  $\mu\text{m}$ , 150  $\times$  4.6 mm I.D.)

Eluent: A) 0.1 % HCOOH in  $\text{H}_2\text{O}$   
 B) 0.1 % HCOOH in  $\text{CH}_3\text{CN}$   
 A/B = 60/40 - 10 min - 0/100, v/v

Flow rate: : 1.0 mL/min


Col. Temp.: : 40  $^\circ\text{C}$

Detection: : UV 214 nm

Injection Vol.: : 10  $\mu\text{L}$

Sample: : ( $\pm$ )-Carboxy-THC 10  $\mu\text{g/mL}$

#### Chemical Structure




( $\pm$ )-11-nor-9-Carboxy- $\Delta^9$ -THC  
(THC-COOH)

Structures are created using Chemistry 4-D Draw which is provided by ChemInnovation Software, Inc.

### What is MonoSpin ?

MonoSpin is a series of spin columns for solid phase extraction (SPE). Owing to the high permeability of monolithic silica disk packed into the spin column, the procedures, such as conditioning, sample loading, washing, and elution can be carried out only by centrifuging the column. It is also the advantage that the elution volume is only 200  $\mu\text{L}$ .

#### An Example of Procedures




## Sample Pretreatment using MonoSpin C18 - CX

MonoSpin C18 - CX has octadecyl group and cation - exchange group on the surface of its silica monolithic support. In this case, acidic buffer was used for sample loading. THC - COOH, which hardly dissociates under acidic condition, could be retained with hydrophobic interaction.


Sample solution	
0.1 % formic acid in water	490 $\mu$ L
Sample matrix	490 $\mu$ L
Standard solution *	10 $\mu$ L
+) Solution of internal standard*	10 $\mu$ L
	1000 $\mu$ L

\* Their concentration is described in the results.



\*\* All the centrifugation was operated at 3,000 xg.

## Evaluation of Recovery by Standard Addition



### Conditions:

<b>System</b>	: LC800 (GL Sciences) API3000 (AB SCIEX)
<b>Column</b>	: InertSustain C18 (2 $\mu$ m, 100 $\times$ 2.1 mm I.D.)
<b>Eluent</b>	: A) 0.1% HCOOH in H <sub>2</sub> O B) 0.1% HCOOH in CH <sub>3</sub> CN A/B = 35/65 - 6 min - 5/95 - 4 min - 5/95 - 0.1 min - 35/65 - 5 min - 35/65, v/v
<b>Flow Rate</b>	: 200 $\mu$ L/min
<b>Col. Temp.</b>	: 40 $^{\circ}$ C
<b>Detection</b>	: MS/MS (ESI, Nega, MRM)
<b>Injection Vol.</b>	: 5 $\mu$ L

\*\*\* Deuterium labeled compound (( $\pm$ )-11-nor-9-Carboxy- $\Delta^9$ -THC-D<sub>9</sub>, THC-COOH-d<sub>9</sub>) was used as internal standard.


### Scan Parameters

Analyte	Precursor Ion ( $m/z$ )	Product Ion ( $m/z$ )
( $\pm$ )-Carboxy-THC	343.0	299.0
( $\pm$ )-Carboxy-THC-d <sub>9</sub>	352.1	308.1


## Calibration Curves

Calibration curve was plotted 4 times by conducting same experiment on different days. All the calibration curve showed good linear response ( $r^2 = 0.9998 \sim 1.000$ ).


1<sup>st</sup> day




2<sup>nd</sup> day



3<sup>rd</sup> day



4<sup>th</sup> day



## Recovery, Reproducibility, and Accuracy

From the chromatograms, recovery, reproducibility, and accuracy were calculated for each day. Almost all results were satisfactory.


Recovery (%)	1 <sup>st</sup> day	2 <sup>nd</sup> day	3 <sup>rd</sup> day	4 <sup>th</sup> day	Inter-day average
0.5 ng/mL	86.8	94.4	87.0	88.5	89.2
1.0 ng/mL	99.0	100.8	101.4	101.7	100.7
10 ng/mL	88.1	91.1	91.6	92.8	90.9

Reproducibility (%)	1 <sup>st</sup> day	2 <sup>nd</sup> day	3 <sup>rd</sup> day	4 <sup>th</sup> day	Inter-day average
0.5 ng/mL	2.67	0.74	2.17	3.20	2.2
1.0 ng/mL	1.00	2.85	2.79	0.57	1.8
10 ng/mL	1.14	1.05	1.79	0.80	1.2

Accuracy (%)	1 <sup>st</sup> day	2 <sup>nd</sup> day	3 <sup>rd</sup> day	4 <sup>th</sup> day	Inter-day average
0.5 ng/mL	101.0	105.3	101.0	101.7	102.2
1.0 ng/mL	99.2	97.1	99.5	98.9	98.7
10 ng/mL	99.9	99.9	99.9	100.1	100.1


## The series of MonoSpin

### MonoSpin C18




Octadecyl group is chemically bonded, and non-polar compounds can be retained because of its hydrophobic interaction. It can be used for extraction or desalting.

### MonoSpin SAX




Trimethylaminopropyl group is bonded. It offers strong anion-exchange and weak hydrophobic interaction. It is suitable for extraction of acidic drugs.

### MonoSpin NH<sub>2</sub>




Aminopropyl group is bonded. It is suitable for extraction of hydrophilic compounds, such as sugar chain.

### MonoSpin PBA




Phenylboronic acid is chemically bonded. Compounds containing *cis*-diol group can be retained with high selectivity.

### MonoSpin SCX




Propylbenzenesulfonic acid is bonded. It offers strong cation-exchange and hydrophobic interaction. It is suitable for extraction of basic drugs.

### MonoSpin TiO




Monolithic silica is coated with titanium dioxide. It is suitable for extraction of phosphate-containing compounds.

### MonoSpin Amide




Amide group is bonded. Wide variety of hydrophilic compounds, such as sugar chain, can be retained and extracted in HILIC mode.

### MonoSpin CBA




Carboxyl group is bonded to silica monolith. It is suitable for extraction of basic compounds through weak cation-exchange interaction.

### MonoSpin C18-CX



Functional groups of both C18 and SCX are bonded. C18-CX is often superior to C18 or SCX in purification and clean-up of basic drugs in serum and urine.

### MonoSpin Trypsin



Trypsin is immobilized on the surface of silica monolithic support. It enables rapid protein-digesting.

(used in this note)

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

#### GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome  
Shinjuku-ku, Tokyo  
163-1130, Japan

Phone: +81-3-5323-6620  
Fax: +81-3-5323-6621  
Email: [world@glsc.co.jp](mailto:world@glsc.co.jp)  
Web: [www.glsciences.com](http://www.glsciences.com)

#### GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255  
Torrance, CA 90503  
USA

Phone: +1-310-265-4424  
Fax: +1-310-265-4425  
Email: [info@glsciencesinc.com](mailto:info@glsciencesinc.com)  
Web: [www.glsciencesinc.com](http://www.glsciencesinc.com)

#### GL Sciences B.V.

Dillenburgstraat 7C  
5652AM, Eindhoven  
The Netherlands

Phone: +31-40-254-9531  
Email: [info@glsciences.eu](mailto:info@glsciences.eu)  
Web: [www.glsciences.eu](http://www.glsciences.eu)

#### GL Sciences (Shanghai) Limited

Tower B, Room 2003  
Far East International Plaza  
No.317 Xianxia Road, Changning District  
Shanghai, China 200051

Phone: +86-21-62782272  
Email: [contact@glsciences.com.cn](mailto:contact@glsciences.com.cn)  
Web: [www.glsciences.com.cn](http://www.glsciences.com.cn)



International Distributors

Visit our Website at [www.glsciences.com/distributors](http://www.glsciences.com/distributors)