




Analysis of iodide ions typically requires separation by ion chromatography followed by detection and quantification using an electrochemical detector (ECD) with a silver electrode, or an inductively coupled plasma mass spectrometer (ICP-MS). However, when using a typical ECD, the silver electrode requires periodic polishing. Both ion chromatography and ICP-MS have the disadvantage of high equipment cost. In this report, an analytical method that combines reversed-phase HPLC and an electrochemical detector with a diamond electrode is shown to be a more convenient method. For separation, iodide ions are successfully retained on an ODS column using a mobile phase with an ion pair reagent. For detection, an electrochemical detector with diamond electrode (that does not require polishing and needs little maintenance) was used. Using this method, it is possible to detect and quantify the iodide ion concentration in seawater.

(C. Aoyama)

**Example: Measurement of standard****HPLC conditions**

<b>Column</b>	: Inertsil ODS-4 (5 μm, 150 x 3.0 mm I.D.)
<b>Eluent</b>	: A) Phosphate buffer * B) CH <sub>3</sub> CN A/B = 65/35, v/v (Mixed by a gradient mixer)
<b>Flow rate</b>	: 0.4 mL/min
<b>Column temperature</b>	: 35 °C
<b>Detected</b>	: ECD 1000 mV vs. Ag/AgCl
<b>Injection volume</b>	: 10 μL


\* Phosphate buffer:  
1.56 g disodium dihydrogen phosphate (dihydrate) 3.58 g  
Disodium hydrogen phosphate (12 hydrate) and 0.32 g  
chloride Hexadecyltrimethylammonium were dissolved in 1  
L of ultrapure water.

**Relationship between potential and peak area****Calibration curve**


## <Example Chromatograms of sea water analysis>

The seawater was filtered through a 0.45  $\mu\text{m}$  membrane filter, which served as a sample solution.

### Seawater sampled at an urban fishing port



### Seawater sampled at a tourist beach



GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

#### GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome  
Shinjuku-ku, Tokyo  
163-1130, Japan

Phone: +81-3-5323-6620  
Fax: +81-3-5323-6621  
Email: [world@glsc.co.jp](mailto:world@glsc.co.jp)  
Web: [www.glsciences.com](http://www.glsciences.com)

#### GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255  
Torrance, CA 90503  
USA

Phone: +1-310-265-4424  
Fax: +1-310-265-4425  
Email: [info@glsciencesinc.com](mailto:info@glsciencesinc.com)  
Web: [www.glsciencesinc.com](http://www.glsciencesinc.com)

#### GL Sciences B.V.

Dillenburgstraat 7C  
5652AM, Eindhoven  
The Netherlands

Phone: +31-40-254-9531  
Email: [info@glsciences.eu](mailto:info@glsciences.eu)  
Web: [www.glsciences.eu](http://www.glsciences.eu)

#### GL Sciences (Shanghai) Limited

Tower B, Room 2003  
Far East International Plaza  
No.317 Xianxia Road, Changning District  
Shanghai, China 200051

Phone: +86-21-62782272  
Email: [contact@glsciences.com.cn](mailto:contact@glsciences.com.cn)  
Web: [www.glsciences.com.cn](http://www.glsciences.com.cn)



**International Distributors**

Visit our Website at [www.glsciences.com/distributors](http://www.glsciences.com/distributors)