


This application describes the analysis of cyclamic acid in foods.

Cyclamic acid (N-cyclohexylsulfamic acid), also called Ticlo, had been used as a sweetener since 1956, but this has been prohibited since 1969, because of reports of suspected carcinogenicity and teratogenicity.

The “test method for cyclamic acid” in the Food Safety Inspection Notification No. 0829009 employs an HPLC method with solid-phase extraction. This report demonstrates that analysis was performed according to this method and excellent data was obtained.

(K.Suzuki, Y.Tanaka)

Example :Measurement of standard






HPLC conditions

Column	: Inertsil ODS-3 (5 μ m, 150 x 4.6 mm I.D.)
Eluent	: A) CH ₃ CN B) H ₂ O A/B = 70/30, v/v (Mixed by a gradient mixer)
Flow rate	: 1.0 mL/min
Column temperature	: 40 °C
Detected	: UV 314 nm
Injection volume	: 10 μ L

Derivatizing reactions

Using hypochlorous acid it is converted to N,N-dichlorocyclohexylamine.



Concentration (mg/L)

Calibration curve

Concentrations indicated before pretreatment.

Example of analysis

Samples pretreated according to the test method described above were used as measurement samples and analyzed by HPLC.

Example of pretreatment

Sample

-Collect 10 g after crushing

Extraction

- Ultrapure water 40 mL
Boiling water (15 min)
Make up to a volume of 100 mL
- Centrifugal separation
Collect-10 mL


Solid phase extraction

- Inertsep Slim C18 (900 mg)
+ Inertsep SAX (500 mg/3 mL)
- Condition each solid phase with 10 mL of methanol
and 10 mL of ultrapure water
- 10 mL sample solution
- Wash with 10 mL of purified water
- Transfer 10mL of hydrochloric acid (1 in 100) to
Inertsep SAX only, and collect.

Derivatized

- To 10 mL of the collected solution add 2 mL of
sulfuric acid (1 in 2) and 5 mL of *n*-hexane,
Add 1 mL of aqueous sodium hypochlorite
solution
- Shake (1min)
- The hexane layer was separated and used as the
sample for measurement.

HPLC-UV



For some samples, the solid-phase cartridge may become clogged.
If this occurs, use a solid-phase cartridge with a large fluid drain area, such as Inertsep C18 (1 g/20 mL).

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly.
We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome
Shinjuku-ku, Tokyo
163-1130, Japan

Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
USA

Phone: +1-310-265-4424
Fax: +1-310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

GL Sciences B.V.

Dillenburgstraat 7C
5652AM, Eindhoven
The Netherlands

Phone: +31-40-254-9531
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (Shanghai) Limited

Tower B, Room 2003
Far East International Plaza
No.317 Xianxia Road, Changning District
Shanghai, China 200051

Phone: +86-21-62782272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn



International Distributors

Visit our Website at www.glsciences.com/distributors