


Simultaneous Analysis of Organic Fluorine Compounds Using LC/MS/MS

As represented by perfluorooctanoic acid (PFOA) and perfluorooctanesulfonic acid (PFOS), organic fluorine compounds are a type of environmental pollutant. These compounds are detected not only in tap water but also in blood of wild animals due to their exceptional stability in the environment. However, their accurate determination has been difficult because these compounds are also eluted from some parts of chromatograph, such as tubing made of PTFE.

In this note, chromatograms obtained with an InertSustain C18 column are shown. Peak shape was excellent for each compound. This analysis was performed on an LC/MS/MS system where the liquid was not in touch with PTFE.





PFBS	-	299/80	4.3
PFHxS	-	399/80	
6.8 PFHxS 13C	-	403/103	
6.8 PFOS	-	499/80	9.9
PFOS 13C	-	503/80	9.9
PFDS	-	599/80	13.1
PFBA	-	213/169	2.2
PFBA 13C	-	217/172	2.2
PFPeA	-	263/219	2.8
PFHxA	-	313/269	
3.7 PFHxA 13C	-	315/270	
3.7 PFHpA	-	363/319	4.8
PFOA	-	413/369	5.9
PFOA 13C	-	417/372	5.9
PFNA	-	463/419	7.3
PFNA	-	468/423	7.3
PFDA	-	513/469	8.8
PFDA	-	515/470	8.7
PFuDA	-	563/519	10.3
PFuDA 13C	-	565/520	10.3
PFdDA	-	613/569	11.9
PFdDA 13C	-	615/570	11.9
PFTTrDA	-	663/619	13.6
PFTeDA	-	713/669	15.1

(2.5 µg/L each)


HPLC Conditions;

System	: LC800
Column	: InertSustain C18 (3 µm, 150 x 2.1 mm I.D.)
Eluent	: A) 10 mM Ammonium acetate B) CH ₃ CN A / B = 60 / 40 – 3 min – 50 / 50 – 22 min – 0 / 100 – 5 min – 0 / 100 (Equilibration for 15 min), v/v
Flow rate	: 0.2 mL/min
Column Temp.	: 40 °C
Detection	: LC/MS/MS (4000Q TRAP®: ESI, Nega, MRM)
Injection Vol.	: 10 µL

* Contact us for more detailed analytical

Chemical structures**PFOA****PFOS**

Structures are created using Chemistry 4-D Draw which is provided by ChemInnovation Software, Inc.



GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan
 22-1 Nishishinjuku 6-chome
 Shinjuku-ku, Tokyo
 163-1130, Japan

Phone: +81-3-5323-6620
 Fax: +81-3-5323-6621
 Email: world@glsc.co.jp
 Web: www.glsciences.com

GL Sciences Inc. USA
 4733 Torrance Blvd. Suite 255
 Torrance, CA 90503
 USA

Phone: +1-310-265-4424
 Fax: +1-310-265-4425
 Email: info@glsciencesinc.com
 Web: www.glsciencesinc.com

GL Sciences B.V.
 Dillenburgstraat 7C
 5652AM, Eindhoven
 The Netherlands

Phone: +31-40-254-9531
 Email: info@glsciences.eu
 Web: www.glsciences.eu

GL Sciences (Shanghai) Limited
 Tower B, Room 2003
 Far East International Plaza
 No.317 Xianxia Road, Changning District
 Shanghai, China 200051

Phone: +86-21-62782272
 Email: contact@glsciences.com.cn
 Web: www.glsciences.com.cn



International Distributors

Visit our Website at www.glsciences.com/distributors