

Reduction of Solvent Consumption in HPLC (2020 #3) –Column Selection

GL Sciences has been suggesting reduction of solvent consumption in HPLC for more than 15 years. This technical note describes how to reduce analysis times and solvent consumption by selection of C18 columns. The example here includes the use of low retentivity C18 columns and changing the column dimensions, which do not require adjustment of the isocratic conditions.


(K. Suzuki)

Low Retentivity C18 Columns

Analysis times and solvent consumption can be reduced only by changing the packing material without any change in conditions or column dimensions. Below is an example of column change from a high retentivity C18 column to InertSustainSwift C18, the retentivity of which is low because of a low carbon load (9%) and a small specific surface area (200 m²/g). If the peaks are separated enough, this technique is recommended as a simple way to reduce the analysis time and solvent consumption.

Analysis of Carotenoids in Seaweed

Column: Competitor C18 (5 μm, 4.6 mm x 150 mm)
Carbon load > 10 %, Specific surface area > 300 m²/g



Conditions


System	: GL7700
Column	: C18 column A (5 μm, 4.6 mm x 150 mm)
Mobile phase	: A) CH ₃ CN B) H ₂ O A/B = 75/25, v/v
Flow rate	: 1.0 mL/min
Column Temp.	: 40 °C
Detection	: VIS 500 nm
Inj. Vol.	: 20 μL
Analyte	: Seaweed extract

1. Fucoxanthin


Note: Please refer to LT030 for the sample pretreatment.

Change to a low retentivity column (same dimensions)

Column: **InertSustainSwift C18** (5 μm, 4.6 mm x 150 mm)
Carbon load: 9%, Specific surface area: 200 m²/g



Comparison of retentivity
among InertSustain and
other brands' C18
columns




Further Reduction of Analysis Time is Possible!


This technique does not require any change in the analytical conditions such as the mobile phase composition. In the example below, a 3-fold reduction of the analysis time is achieved by simply changing the particle diameter and column dimensions. This reduction of the analysis time also decreases the solvent consumption.

3-fold Reduction of Analysis Time and Solvent Consumption!

Column : InertSustainSwift C18
5 μm , 4.6 mm x 150 mm
 Flow rate : 1.0 mL/min



Column : InertSustainSwift C18
3 μm , 3.0 mm x 100 mm
 Flow rate : 1.0 mL/min




Other Features of InertSustainSwift C18

· Fast equilibration!


· For small and mid-sized (up to $M_w \approx 20,000$) molecules!

Comparison of Equilibration



Please contact us for more detailed technical information.

Analysis of 3 oligonucleic acids



Please refer to LB514 for the analytical conditions.

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome
 Shinjuku-ku, Tokyo
 163-1130, Japan

Phone: +81-3-5323-6620
 Fax: +81-3-5323-6621
 Email: world@glsc.co.jp
 Web: www.glsciences.com

GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255
 Torrance, CA 90503
 USA

Phone: +1-310-265-4424
 Fax: +1-310-265-4425
 Email: info@glsciencesinc.com
 Web: www.glsciencesinc.com

GL Sciences B.V.

Dillenburgstraat 7C
 5652AM, Eindhoven
 The Netherlands

Phone: +31-40-254-9531
 Email: info@glsciences.eu
 Web: www.glsciences.eu

GL Sciences (Shanghai) Limited

Tower B, Room 2003
 Far East International Plaza
 No.317 Xianxia Road, Changning District
 Shanghai, China 200051

Phone: +86-21-62782272
 Email: contact@glsciences.com.cn
 Web: www.glsciences.com.cn



International Distributors

Visit our Website at www.glsciences.com/distributors