



Analysis of glucosamine using a differential index detector

Glucosamine is an amino sugar in which one of the hydroxy groups of glucose is replaced with an amino group, and is often used in health foods. At the 10th Edition Food Additive Official Formulation Study Group (3rd), a draft ingredient standard using liquid chromatography as a quantification method for glucosamine was decided. This time, we would like to introduce the analysis of glucosamine that complies with the draft component standard using the Hitachi HPLC system Chromaster. (K. Tanaka)

Chemical Structures




Standard solution measurement example



Conditions

- Column** : InertSustain NH2
(5 μm, 250 x 4.6 mm I.D.)
- Eluent** : A) CH₃CN
B) H₂O
A/B = 3/1, v/v
- Flow Rate** : 0.7 mL/min *
* Adjust the flow rate so that the retention time of glucosamine is about 12 minutes.
- Col. Temp.** : 40 °C
- Detection** : RI
- Injection Vol.** : 10 μL

Glucosamine calibration curve (0~10 mg/mL)



Food additives, such as preservatives, sweeteners, colorants, and flavors, are used in the manufacturing process of foods or for the purpose of processing and preserving foods. The Ministry of Health, Labor and Welfare has received an evaluation by the Food Safety Commission regarding the safety of food additives, and only when there is no risk of damaging human health, the Ministry of Health, Labor and Welfare approves the use after establishing the standard of ingredients and the standard of use. I am. The official food additive standard stipulates standards for ingredients of food additives, manufacturing standards, and methods for ensuring quality. (From the Ministry of Health, Labor and Welfare website)

Product used

●Analytical column

InertSustain NH2 with 100% CH₃CN

5 μm, 250 x 4.6 mm I.D.

Cat.No. 5020-89951

☆ This is a product (usually hexane / ethanol type) in which the encapsulation solvent at the time of shipment is changed to 100% acetonitrile. Recommended if you want to use it in HILIC mode.



GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome
Shinjuku-ku, Tokyo
163-1130, Japan

Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
USA

Phone: +1-310-265-4424
Fax: +1-310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

GL Sciences B.V.

Dillenburgstraat 7C
5652AM, Eindhoven
The Netherlands

Phone: +31-40-254-9531
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (Shanghai) Limited

Tower B, Room 2003
Far East International Plaza
No.317 Xianxia Road, Changning District
Shanghai, China 200051

Phone: +86-21-62782272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn



International Distributors

Visit our Website at www.glsciences.com/distributors