


Analysis of sugars with an InertSphere Sugar-2 Column -Using Ligand-Exchange and Size-Exclusion Modes-

InertSphere Sugar-2 is a sugar analysis column that uses calcium (Ca^{2+} type) as the metal counter ion in the packing material. The separation mainly occurs due to a size-exclusion mechanism and sugars elutes in the order of molecular weight. It also acts in a ligand exchange mode by utilizing the difference in retention between the metal counter ions and hydroxyl groups in the sugars. The bond strength of the complexation depends on the counter ion and the type of sugar.

One of the advantages of this column is that a 100% aqueous eluent is used, eliminating the need for sample preparation.


The InertSphere Sugar-2 column was used to analyze sugars in vegetable juices and yogurts. As a comparison, similar analyses were performed using an amino column (NH_2), separation of difficult components was achieved.. (Y. Yui)

Example: Measurement of standards



Typically NH_2 columns cannot separate the overlapping peaks of Glucose, Galactose, and Sorbitol,



Using an InertSphere Sugar-2 column, all five components could be separated.





Calibration curve (0 mg/mL - 100 mg/mL)

HPLC-RI Conditions
Column : InertSphere Sugar-2
 (9 μm , 300 x 7.8 mm I.D.)
Eluent : H_2O
Flow Rate : 0.5 mL/min
Col. Temp. : 85 °C
Detection : RI
Injection Vol. : 10 μL
Sample : Standard

Examples of pretreatment and analysis of sugars in vegetable juice.



Examples of pretreatment and analysis of sugars from plain yogurt.



HPLC Conditions

Guard column : InertSphere Sugar-2 Guard
(9 µm, 50 x 6.0 mm I.D.)
Column : InertSphere Sugar-2
(9 µm, 300 x 7.8 mm I.D.)
Eluent : H₂O

Flow Rate : 0.5 mL/min
Col. Temp. : 85 °C
Detection : RI
Injection Vol. : 10 µL

Retention time list (reference value)

Natural sugar	
Glucose	12.66
Galactose	13.98
Xylose	13.90
Mannose	14.37
Arabinose	15.69
Fructose	15.45
Ribose	24.10
Disaccharide	
Trehalose	10.50
Sucrose	10.52
Maltose	10.61
Lactose	10.90
Palatinose	10.81
Isomaltose	10.48
Sugar alcohol	
Dulcitol	22.55
Sorbitol	23.55
Myo-inositol	15.75
Xylitol	23.39
Arabitol	20.12
Maltol	65.70
Maltitol	13.99
Lactitol	13.81
Mannitol	19.63
Rare seven monosaccharides	
Glucoheptose	17.24
Rare six monosaccharides	
Sorbose	14.12
Growth	17.58
Tagatose	18.38
Talose	22.09
Fucose	15.66
Rhamnose	14.54
Allulose	18.32
Rare pentasaccharide	
Lyxose	16.39
Rare tetramonosaccharide	
Threose	15.51
Monosaccharide	
N-acetyl-D-glucosamine	13.42
Oligosaccharides	
1-Kestose	9.42
1F-fructofuranosyl nystose	8.63
Nystose	8.90
Maltotetraose	9.18
Maltotriose	9.69
Maltohexaose	8.64
Maltoheptaose	8.50
Maltopentaose	8.59
Meretitose	9.52
Raffinose	9.62
Isomaltotriose	9.46
Amino sugar	
Glucosamine	8.44
Galactosamine	8.45
Artificial sweetener	
Aspartame	14.24
Acesulfam K	8.57
Advantame	20.24
Saccharin	20.16

HPLC Conditions

Column	: InertSphere Sugar-2 (9 µm, 300 x 7.8 mm I.D.)
Eluent	: H ₂ O
Flow Rate	: 0.5 mL/min
Col. Temp.	: 85 °C
Detection	: RI
Injection Vol.	: 10 µL
Sample	: Standard

Column

**Analytical columns: InertSphere Sugar-2 9 µm, 300 x 7.8 mm I.D.
Cat.No. 5020-11000**

**Guard columns: InertSphere Sugar-2 Guard 9 µm, 50 x 6.0 mm I. D.
Cat.No. 5020-10999**

- Base Material : Styrene-divinylbenzene-based polymer
- Particle size : 9 µm
- Functional Group : Sulfonic acid group
- Counter-ion : Ca²⁺
- Degree of linking : 8 %
- USP code : L19



GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome
Shinjuku-ku, Tokyo
163-1130, Japan

Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
USA

Phone: +1-310-265-4424
Fax: +1-310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

GL Sciences B.V.

Dillenburgstraat 7C
5652AM, Eindhoven
The Netherlands

Phone: +31-40-254-9531
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (Shanghai) Limited

Tower B, Room 2003
Far East International Plaza
No.317 Xianxia Road, Changning District
Shanghai, China 200051

Phone: +86-21-62782272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn



International Distributors

Visit our Website at www.glsciences.com/distributors