

InertSearch™ for LC

Inertsil® Applications

Simultaneous analysis of bile acids using LC/MS/MS

Data No. LL030-0000

The chromatograms were provided by Masamitsu Maekawa and Nariyasu Mano, Department of Pharmaceutical Sciences, Tohoku University Hospital, 1-1 Seiryō-machi, Aoba-ku, Sendai, Miyagi, 980-8574, Japan

Conditions

Column : InertSustain C18
(2 μm, 150 x 2.1 mm I.D.)
Column Cat. No. : 5020-14355
Eluent : A) 0.1% HCOOH in H₂O
B) CH₃CN/CH₃OH/HCOOH = 50/50/0.1, v/v/v

Time(min)	A (vol%)	B (vol%)
0	60	40
10	45	55
20	0	100

Flow rate : 0.4 mL/min
Col. Temp. : 40 °C
Detection : MS/MS (ESI, Negative, MRM)
Injection Vol. : 5 μL
Sample : Standard solution

Analyte	Precursor ion (m/z)	Product ion (m/z)
1. Taurine-conjugated α-muricholic acid (TaMCA)	514.2	80.0
2. Taurine-conjugated ω-muricholic acid (TωMCA)	514.2	80.0
3. Taurine-conjugated β-muricholic acid (TbMCA)	514.2	80.0
4. Taurine-conjugated γ-muricholic acid (THCA)	514.2	80.0
5. Taurine-conjugated ursodeoxycholic acid (TUDCA)	498.2	80.0
6. Taurine-conjugated hyodeoxycholic acid (THDCA)	498.2	80.0
7. Glycine-conjugated α-muricholic acid (GaMCA)	464.2	74.0
8. Glycine-conjugated β-muricholic acid (GbMCA)	464.2	74.0
9. Taurine-conjugated cholic acid (TCA)	514.2	80.0
10. Glycine-conjugated γ-muricholic acid (GHCA)	464.2	74.0
11. Glycine-conjugated ursodeoxycholic acid (GUDCA)	448.2	74.0
12. ω-Muricholic acid (ωMCA)	453.2	407.2
13. Glycine-conjugated hyodeoxycholic acid (GHDCA)	448.2	74.0
14. α-Muricholic acid (αMCA)	453.2	407.2
15. Taurine-conjugated chenodeoxycholic acid (TCDCA)	498.2	80.0
16. β-Muricholic acid (βMCA)	453.2	407.2
17. Glycine-conjugated cholic acid (GCA)	464.2	74.0
18. Taurine-conjugated deoxycholic acid (TDCA)	498.2	80.0
19. γ-Muricholic acid (HCA)	453.2	407.2
20. Ursodeoxycholic acid (UDCA)	437.2	391.2
21. Hyodeoxycholic acid (HDCA)	437.2	391.2
22. Taurine-conjugated lithocholic acid (TLCA)	482.2	80.0
23. Cholic acid (CA)	453.2	407.2
24. Glycine-conjugated chenodeoxycholic acid (GCDCA)	448.2	74.0
25. Glycine-conjugated deoxycholic acid (GDCA)	448.2	74.0
26. Chenodeoxycholic acid (CDCA)	437.2	391.2
27. Glycine-conjugated lithocholic acid (GLCA)	432.2	74.0
28. Deoxycholic acid (DCA)	437.2	391.2
29. Lithocholic acid (LCA)	421.3	391.2

InertSearch™ for LC

Inertsil® Applications

Simultaneous analysis of bile acids using LC/MS/MS

